Instructions for Use

Bulgarian (BG)

Croatian (HR)

Czech (CS)

Danish (DA)

English-UK (GB)

Estonian (ET)

French (FR)

German (DE)

Greek (EL)

Hungarian (HU)

Italian (IT)

Lithuanian (LT)

Macedonian (MK)

Norwegian (NO)

Polish (PL)

Portuguese (PT)

Romanian (RO)

Russian (RU)

Slovak (SK)

Spanish (ES)

Swedish (SV)

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

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
Nederlandene

Exact Diagnostics LLC
3400 Camp Bowie Blvd CBH-214
Fort Worth Texas 76107
USA

EDX Babesia Run Control Summary	
Catalogue number	BABRC
Analyte	<i>Babesia microti</i>
Matrix	Synthetic Matrix
Preservative	ProClin® 300
Storage	-20 °C or below
Fill volume	1.0 mL
Number of vials per kit	5 vials
Number of uses	Multiple use
Shelf life	18 months from date of manufacture
Precautions	Biological risks
Regulatory status	For in vitro diagnostic use

INTENDED USE:

EDX Babesia Run Control is an external quality control intended to be used with *Babesia microti* assays to monitor the performance of that given assay. The control is used to monitor the presence of *Babesia microti* DNA. Routine use of EDX Babesia Run Control allows laboratories to evaluate day-to-day and lot-to-lot variation of their molecular assay and test for operator proficiency.

PRODUCT DESCRIPTION:

EDX Babesia Run Control contains whole, intact organisms. The control is intended to test the entire process of a molecular assay including extraction, detection, and amplification.

EDX Babesia Run Control is internally value assigned using Bio-Rad Digital Droplet PCR (ddPCR). EDX Babesia Run Control is formulated in Synthetic Matrix that contains preservatives.

PROCEDURE:

EDX Babesia Run Control requires an extraction step. The control must be treated in a similar manner to other tested samples. This allows the operator to assess their extraction technique.

EDX Babesia Run Control is an unassayed control and does not have assigned values. The control is not a substitute for the manufacturer’s kit controls provided with the assay. EDX Babesia Run Control is to be tested according to the assay manufacturer’s or testing laboratory’s instructions and recommendations. Expected results using EDX Babesia Run Control must be established by the end user for their specific assay.

STORAGE AND HANDLING:

EDX Babesia Run Control should be stored at -20 °C or below.

Thaw EDX Babesia Run Control at room temperature and vortex briefly prior to use.

Once thawed, EDX Babesia Run Control is stable for 24 hours when stored at 2-8 °C. Product can be frozen and reused a maximum of four (4) times. After the fourth freeze-thaw, any remaining materials must be disposed of. Do not dilute.

Do not use EDX Babesia Run Control beyond the expiration date.

LIMITATIONS:

For in vitro diagnostic use.

WARNINGS AND PRECAUTIONS:

EDX Babesia Run Control should be considered biohazardous. Universal precautions and proper disposal should be practised¹. Do not pipette by mouth. Do not smoke, eat, or drink in areas where specimens are handled.

Discard product if packaging is damaged or leaking. Disinfect liquids, materials, or spills with a 0.5% sodium hypochlorite solution. Dispose of all materials and liquids used in the procedure as if they contain pathogenic agents.

REFERENCES:

¹ Siegel JD, Rhinehart E, Jackson M, Chiarello L, and the Healthcare Infection Control Practices Advisory Committee, 2007
Guideline for Isolation Precautions: Preventing transmission

BIO-RAD and DROPLET DIGITAL PCR are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions.

SYMBOL REFERENCES:

Catalogue number

Lot number

For in vitro diagnostic use

Biological risks

Expiration date

Upper limit of temperature

Manufacturer

Caution

Consult instructions for use


Positive Control

European Mark of Conformity

Authorised Representative

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Exact Diagnostics LLC
3400 Camp Bowie Blvd CBH-214
Fort Worth, Texas 76107
United States

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BIO-RAD och DROPLET DIGITAL PCR är varumärken som tillhör Bio-Rad Laboratories, Inc. i vissa länder.