

EDX CMV LOD Summary	
Catalog Number	CMVLOD
Analyte	Inactivated CMV
Matrix	True EDTA Plasma
Preservative	ProClin® 300
Storage	-70°C or Below
Fill Volume	1.0mL
Number of Vials Per Kit	2 x 10 Vials
Number of Uses	Single Use
Shelf Life	24 Months from Date of Manufacture
Traceability	1 <sup>st</sup> WHO International Standard
Precautions	Biological Risks
Regulatory Status	For In Vitro Diagnostic Use

**INTENDED USE:**

The EDX CMV LOD is calibrated against the 1<sup>st</sup> WHO International Standard for Cytomegalovirus (09/162)<sup>1</sup>. EDX CMV LOD is used to determine the analytical sensitivity of the assay for Cytomegalovirus (CMV) DNA. EDX CMV LOD allows laboratories to evaluate their molecular assay and test for operator proficiency.

**PRODUCT DESCRIPTION:**

The EDX CMV LOD contains whole, intact virus. EDX CMV LOD is developed using heat inactivated CMV Merlin strain. The control is intended to test the entire process of a molecular assay including extraction, detection, and amplification.

The EDX CMV LOD is formulated in True EDTA Plasma that is non-reactive for HBV DNA, HCV RNA, HIV-1 RNA, CMV DNA, HBsAg and antibodies for HIV-1/2, and HCV.

**PROCEDURE:**

The EDX CMV LOD requires an extraction step. The control must be treated in a similar manner to other tested samples. This allows the operator to assess their extraction technique.

The EDX CMV LOD is an unassayed control and does not have assigned values. The control is not a substitute for the manufacturer's kit controls provided with the assay. EDX CMV LOD is to be tested according to the assay manufacturer's or testing laboratory's instructions and recommendations. Expected results utilizing EDX CMV LOD must be established by the end user for their specific assay.

**STORAGE AND HANDLING:**

The EDX CMV LOD should be stored at -70°C or below.

Thaw EDX CMV LOD at room temperature and vortex briefly prior to use.

Once thawed, EDX CMV LOD is stable for 24 hours when stored at 2-8°C. Any remaining materials must be disposed of after one use. Do not reuse. Do not dilute.

Do not use EDX CMV LOD beyond the expiration date.

**LIMITATIONS:**

For In Vitro Diagnostic Use.

**WARNINGS AND PRECAUTIONS:**

The EDX CMV LOD contains heat inactivated virus but should be considered biohazardous. Universal precautions and proper disposal should be practiced<sup>2</sup>. Do not pipette by mouth. Do not smoke, eat, or drink in areas where specimens are handled.

Discard product if packaging is damaged or leaking. Disinfect liquids, materials or spills with a 0.5% sodium hypochlorite solution. Dispose of all materials and liquids used in the procedure as if they contain pathogenic agents.

**REFERENCES:**

<sup>1</sup> The 1<sup>st</sup> WHO International Standard for Cytomegalovirus (09/162) IFU. <http://www.nibsc.org/documents/ifu/09-162.pdf>

<sup>2</sup> Siegel JD, Rhinehart E, Jackson M, Chiarello L, and the Healthcare Infection Control Practices Advisory Committee, 2007 Guideline for Isolation Precautions: Preventing Transmission

**SYMBOL REFERENCES:**



Catalog Number



Lot Number



For In Vitro Diagnostic Use



Biological Risks



Expiration Date



Do not re-use



Upper Limit of Temperature



Manufacturer



Caution



Consult Instructions for Use



Positive Control



Exact Diagnostics LLC  
3400 Camp Bowie Blvd CBH-214  
Fort Worth TX 76107  
United States