

PRODUCT INFORMATION

Product SKU	HSVP200
Product Name	HSV Verification Panel
Kit Lot Number	19112143
Source Material 1	Inactivated HSV-1
Source Material 2	Inactivated HSV-2
Date of Expiration	2021-05-31
Base Matrix	Multiprep Specimen Diluent
Preservative	Sodium Azide
Fill Volume	1.0 mL
Storage Temperature	-20°C or Below
Regulatory Status	CE-IVD
Testing Platform	cobas® HSV 1 and 2 Test

This is an unassayed control and does not have assigned analyte values. The end user is responsible for determining the value of the control. This quality control material should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols.

TESTING DATA

Lot Number	Product Description	Results
19112134	HSV-1 3x LOD	Pass
19112135	HSV-1 5x LOD	Pass
19112136	HSV-1 10x LOD	Pass
19112137	HSV-2 3x LOD	Pass
19112138	HSV-2 5x LOD	Pass
19112139	HSV-2 10x LOD	Pass
19112140	HSV-1/HSV-2 3x LOD	Pass
19112141	HSV-1/HSV-2 5x LOD	Pass
19112142	HSV-1/HSV-2 10x LOD	Pass

Bioburden Testing Results	
Acceptance Criteria: Less than 10 CFU/mL	Pass

APPROVAL

Quality Assurance	Date
<i>Taylor Taylor</i>	03Jan20